

**Low Foam Detergent**

Product information Data Sheet

**PRODUCT DESCRIPTION:**

*Low Foam Detergent* is a highly concentrated neutral detergent and is safe for use on all types of instrument and equipment.

*Low Foam Detergent* is recommended for instruments that have already been pre-cleaned.

*Low Foam Detergent* has a built-in lubrication agents which prolongs instruments life.



**RECOMMENED USES:**

- Medical/surgical instruments
- Stainless steel and aluminum items
- Respiratory Therapy Equipment
- Evacuation Systems
- Tungsten and carbide copper items

**ADVANTAGES:**

*Low Foam Detergent*

- Safe and effective for stainless steel, aluminum, tungsten, carbide copper, and brass
- Improves cleaning and rinses freely
- Formulated for all brands of automatic washer and sterilizers
- Built-in lubrication agents to prolong instrument life.

**TECHNICAL DATA:**

- pH neutral
- Non-corrosive
- No Phosphates, Chlorides or Hydroxides
- Dilution 1:128 (1 ounce per gallon)

**PACKAGING:**

- 1 Gallon
- 5 Gallons
- 15 Gallons
- 30 Gallons

MANUFACTURES DATA AND SAFETY REQUIREMENTS			
Specific Gravity	1.01	Safety Equipment and Clothing	Goggles, Gloves, and Coveralls
Boiling Temperature	108 <sup>0</sup> C	Disposal Guidelines	See MSDS
Storage Temperature	Room Temperature	Flash Point	None

Manufactured by:  
Enzyme Solutions, Incorporated  
7601 Honeywell Drive  
Fort Wayne, Indiana 46825  
260.497.0851 Telephone  
260.497.0512 Fax  
[www.enzymesolutions.com](http://www.enzymesolutions.com)