

# Neutralizing Rinse

## Product information Data Sheet

**PRODUCT DESCRIPTION:**  
*Neutralizing Rinse* is a mild acidic solution designed to neutralize the alkaline residue that detergents leave on glassware, surgical instruments, rubber, metal, and plastic items.  
*Neutralizing Rinse* leaves no spots and keeps instruments and equipment looking like new.  
*Neutralizing Rinse* ensures complete neutralization of residual alkaline detergent carried over from the cleaning cycle in automatic washer/decontaminators.



**RECOMMENDED USES:**

- Medical/surgical instruments
- Glassware
- Rubber, metal, and plastic items
- Respiratory Therapy equipment

**ADVANTAGES:**

*Neutralizing Rinse*

- Acidic rinse promotes natural passivation process for additional instruments protection
- The anti-corrosive agents inhibits corrosion which prolongs the instruments use life
- Chelating additives counter act hard water and help protect surgical instruments by preventing spotting and staining
- Super concentrated to save money

**TECHNICAL DATA:**

- pH 1.7
- Non-flammable
- Dilution 1-3 ounces per rinse cycle

**PACKAGING:**

- 5 Gallons
- 15 Gallons

MANUFACTURES DATA AND SAFETY REQUIREMENTS			
Specific Gravity	1.17	Safety Equipment and Clothing	Goggles, Gloves, and Coveralls
pH	1.7	Disposal Guidelines	See MSDS
Storage Temperature	Room Temperature	Flash Point	None

Manufactured by:  
 Enzyme Solutions, Incorporated  
 7601 Honeywell Drive  
 Fort Wayne, Indiana 46825  
 260.497.0851 Telephone  
 260.497.0512 Fax  
[www.enzymesolutions.com](http://www.enzymesolutions.com)