

Low Foam Detergent

Low Foam Detergent is a highly concentrated neutral detergent and is safe for use on all types of instrument and equipment.

Low Foam Detergent is recommended for instruments that have already been pre-cleaned.

Low Foam Detergent has a built-in lubrication agents which prolongs instruments life.



RECOMMENDED USES:

- Medical/surgical instruments
- Stainless steel and aluminum items
- Respiratory Therapy Equipment
- Evacuation Systems
- Tungsten and carbide copper items

ADVANTAGES:

- Safe and effective for stainless steel, aluminum, tungsten, carbide copper, and brass
- Improves cleaning and rinses freely
- Formulated for all brands of automatic washer and sterilizers
- Built-in lubrication agents to prolong instrument life.

TECHNICAL DATA:

- pH neutral
- Non-corrosive
- No Phosphates, Chlorides or Hydroxides
- Dilution 1/8 to 1 ounce per gallon

PACKAGING:

- 4x1 Gallon
- 2x2.5 Gallons
- 5 Gallons
- 15 Gallons
- 30 Gallons

ENZYME SOLUTIONS®
INCORPORATED

2105 Forrest Park Drive
Garrett, IN 46738
800-523-1323 Phone • 260.497.0512 Fax
www.enzymesolutions.com