

Low Foam Detergent **ULTRA CONCENTRATE**

Low Foam Detergent ULTRA CONCENTRATE is a highly concentrated neutral detergent and is safe for use on all types of instrument and equipment.

Low Foam Detergent ULTRA CONCENTRATE is recommended for instruments that have already been pre-cleaned.

Low Foam Detergent ULTRA CONCENTRATE has built-in lubrication agents which prolong instrument life.



RECOMMENDED USES:

- Medical/surgical instruments
- Stainless steel and aluminum items
- Respiratory Therapy Equipment
- Evacuation Systems
- Tungsten and carbide copper items

ADVANTAGES:

- Safe and effective for stainless steel, aluminum, tungsten, carbide copper, and brass
- Improves cleaning and rinses freely
- Formulated for all brands of automatic washers and sterilizers
- Built-in lubrication agents to prolong instrument life

TECHNICAL DATA:

- pH neutral
- Non-corrosive
- No Phosphates, Chlorides or Hydroxides
- Dilution 1/20 to 1/10 ounce per gallon

PACKAGING:

- 4x1 Gallon
- 2x2.5 Gallons

ENZYMESOLUTIONS®
INCORPORATED

2105 Forrest Park Drive
Garrett, IN 46738
800-523-1323 Phone • 260-497-0512 Fax
www.enzymesolutions.com